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
Decontaminates equipment and PPEs contaminated by chemical products.

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PREVOR
ANTICIPATE AND SAVE
Toxicology Laboratory & Chemical Risk Management

DEVICE IN CASE OF SPILLS OF CHEMICAL PRODUCTS

The C.M.I.C. 78 (Yvelines Chemical Intervention Mobile Unit) has equipped itself with the TRIVOREX® neutralizing absorbent (decontamination of chemical spills) in its operational capabilities since a few years.
PREVOR organized a chemical risks training session with the C.M.I.C. 78 incorporating the problematic of chemical spills. Fifteen firemen "RCH3" and "RCH4" qualified participated. We reminded equipment and visualized the procedure (colored indicator of integrated neutralization) dealing with liquid spills with or without gas release. This also allowed us to exchange experience feedbacks from industrials and interventions like:

- ✓ Spillage of a battery (sulfuric acid) on a barge treated by neutralization and simultaneous absorption of the TRIVOREX® neutralizing absorbent.
- ✓ Incompatible mix of hydrochloric acid (HCl) and bleach (sodium hypochlorite NaClO) in a water treatment technical premise of a swimming pool with release of toxic chlorinated gas. The use of the TRIVOREX® neutralizing absorbent brought back the liquid to a neutral pH with a stop of the toxic vapors release and the security within the swimming pool.

The Commandant Aymeric ARNOULD, Yvelines Departmental Chemical Risks Technical adviser:

"TRIVOREX® is polyvalent, that can be used and is adapted to most of our interventions we can be fronted to.

Its practical aspect presents an advantage linked to the simplicity of implementation without necessity of neutralization calculations.

From this practical aspect follows a security aspect equally important during an intervention."

Editor : Goulven TILLY
Safety Environment Engineer
PREVOR

Testimony : Commandant ARNOULD
Yvelines Departmental Chemical Risks
Technical adviser – C.M.I.C. 78

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PREVOR ENVIRONMENT
HANDLING AND RESTORATION
Toxicology laboratory and chemical risk management

Courtesy translation of the SDIS 42 (French fireman brigade) protocol of use for TRIVOREX® neutralizing absorbent.

FOS CBRN
Technical sheet RCH
N°033

« Application of Trivorex® neutralizing powder »


Redactor: LTN REYMOND Date: 26/01/2015

Definition-description TRIVOREX® absorbent is a powder to absorb, solidify and neutralize spills of liquid chemical product.

Indications: The powder is packed inside a 10 kg bucket.

The powder color changes by itself depending on the pH of the product to be treated.

Absorbent with integrated pH indicator

When the neutralization is achieved, the mix turns into a gel that will have to be collected and disposed in accordance with the specific waste treatment process borne by the polluter.

Contraindications: This powder is not used on water.

Application

⚡ **Can be used on all types of :**

- Acid
- Bases
- Hydrocarbons

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Principles of use :

It is mandatory to wear PPEs defined by the CBRN officer.

Maintenance-Renewal Make a replacement request through the intranet.

Miscellaneous This product can also be used on organic waste

Courtesy translation of the SDIS 42 (French fireman brigade) decontamination protocol of chemical protection suits.

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Definition- description	Decontamination solution for suits to use during interventions of chemical nature (acids and bases).
------------------------------------	--

Indications:	<p>The solution is packed inside a 10 liters can.</p>  <p>The solution must be used pure. Depending on the pH of the contaminant, the neutralization is more or less quick (around one minute).</p>
---------------------	---

contraindication:	Keep product from freezing.
--------------------------	-----------------------------

Application	<p>+ Usable on all types of:</p> <ul style="list-style-type: none"> • Acids • Bases <p>+ Principles of use:</p> <ul style="list-style-type: none"> - This action of decontamination will be made inside the decontamination chamber as soon as the pairs partner return.  <p>The solution will be pulverized directly on the chemical suit starting with the head.</p> <ul style="list-style-type: none"> - Start with the head of the wearer. - Pulverize on the entire chemical suit surface. - Wait at least one minute.
--------------------	--

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	<ul style="list-style-type: none"> - Rinse the chemical suit with water. - In case of doubt, rub a pH paper band on the chemical suit litigious areas.  <ul style="list-style-type: none"> - Proceed to the classical undressing of the wearer.
--	---

Maintenance- Renewal	Make a replacement request through the intranet.
---------------------------------	--

Various:	- This product can also be used on organic waste.
-----------------	---

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OBI Andreas Kapeller, MSc

Linz, am 19.12.2008

andreas.kapeller@btf-chemiepark.com

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